1

5

2

6


3

7

4

8

Deutsch

Collamed® long

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Collamed® long

Intended purpose

Collamed is a knee support brace which functions based on the 4-point principle using limitation of extension / flexion.

Indications

All indications that require reinforcement of the physiological stabilisation of the knee joint and/or stress relief of the joint, such as:

- Following injury to the collateral ligament (post-traumatic/postoperative)
- Following injury to the cruciate ligament (post-traumatic/postoperative)
- For knee instability
- Following meniscal injury
- For knee osteoarthritis

Contraindications

None known at present.

Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in

conjunction with sweating) or due to the composition of the material.

Intended patient groups

Healthcare professionals should provide care to the adults and children, applying the available information on the measurements/sizes and necessary functions/indications and in line with the information provided by the manufacturer, and acting under their own responsibility.

Setting the flexion / extension limits Only to be carried out by an orthopedic technician.

Extension limit wedge, EXT (limits extension): 0°, 10°, 20°, 30°

Flexion limit wedge, FLEX (limits flexion): 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°

1. Your knee brace is equipped with 0° extension wedges and 90° flexion wedges.
2. To change the wedges, open the top and bottom sliding lock. Now, the wedges can simply be taken out.
3. The replacement wedges are included in the set of wedges labelled with R - FLEX, L - FLEX or R - EXT, L - EXT and the respective degree settings.
4. Firstly, insert the new flexion wedge and then the extension wedge into the joint. Afterwards, close the two sliders to secure the wedges. Repeat the same process at the other joint.
WARNING: Extension and flexion wedges must always be inserted together.
5. Move the brace joints to check that the wedges are positioned securely.

Important notes

Only change the extension and flexion limits on instruction from your doctor or orthotist.

In order to prevent hyperextension,

always insert extension stops or 0° stops. The stops define the exact excursion of the orthosis. In difficult soft-tissue conditions, e.g. where the soft-tissues are under compression, it may be necessary to insert a larger stop to limit the joint movement to the desired excursion. In addition, make sure that the orthosis hinges are correctly positioned: the middles of the hinges should be in one line with the upper border of the patella. The hinge must lie behind the midline of the leg.

Fitting instructions

- Open all straps and the two front hook and loop fasteners. Release the straps from one of the two hook and loop points.
- Put your leg into a 45-degree angle flexion position. Wrap the brace closely around the leg from behind and close the two hook and loop fasteners (Image 1).
- Make sure that the kneecap is positioned centrally in the space for the patella. Position the brace in such a way that the centres of the hinges are level with the upper edge of the kneecap and just behind the middle of the side of the leg. (Image 2).
- Now, fasten the bottom calf strap (Image 3) and then the top thigh strap (Image 4).
- Close the top calf strap (Image 5) and the bottom thigh strap (Image 6).
- If possible, walk a few steps in order to check whether the brace is in the correct position (Image 4). Please adjust the position as required (Image 7).

Care instructions

Please close the hook and loop fasteners before washing and remove the orthoses. Soap residues, lotions and ointments can cause skin irritation and material wear.

- Wash the product by hand, preferably using medi clean washing agent.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

Storage instructions

Keep the brace in a dry place and do not expose to direct sunlight.

Material composition

Avional, polyamide, elasthane

Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

Disposal

The product can be disposed of in the domestic waste.

Your medi team

wishes you a quick return to full fitness

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR). The traceability of this product is assured via a UDI code [UDI](#).

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