SecuTec[®] Genu Flex

STAYS IN PLACE DURING MOVEMENT

KNFF TRFATMENT

BAUERFEIND YOUR KNEE THERAPY PARTNER

All-round treatment for cruciate ligament, collateral ligament, and meniscus injuries



SofTec[®] Genu







For a complete list of all indications, please see the corresponding instructions for use.

SecuTec[®] Genu Flex WITH SPECIAL SIZING SYSTEM Maximum size coverage with less versions.

SIZE	1	2	3	4
Circumference top (cm)	37-45	45-53	53-61	61 - 69
Circumference bottom (cm)	28-36	32-40	36-44	40-48

ITEM NUMBER 121 4677 1050 +0°

+ = Design (6) right (7) left, ° = Size (1 - 5)

INDICATIONS

- Ruptured anterior / posterior cruciate ligament (ACL / PCL)
 To restrict the knee joint's range of motion
- After ligament surgery / ligament reconstruction
- Severe and / or complex instability
- Collateral ligament injuries

- (e.g. after meniscus surgery)
- Non-surgical treatment or long-term rehabilitation
- Functional instabilities

All information in this brochure is up to date at the time of printing. This brochure – depending on its version – gives an initial insight into this area of expertise, and suggests selected medical diagnosis and treatment options. Even though the details are scientifically substantiated, we do not claim that all information is complete or that all schools of thought are represented. All contents in this brochure or its online version, including extracts, in particular those about diagnostics and treatment options (texts, images, graphs, etc.), are understood as non-binding support for the treatment that must be chosen by a physician or specialist. This brochure is no substitute for the necessary personal medical history and its overall assessment by a physician / a professionally trained and qualified specialist, considering ongoing treatment, if applicable. The contents in this brochure are exclusively designed for medical specialist staff, not for anyone unrelated to the subject, in particular not for patients. If in doubt, please seek additional advice from a specialist. We are also happy to help.

For long-term treatment of complex instabilities.*

VenoTrain[®] soft S

b:joynz knee immobilizer 0/20 For immobilization of the knee.

SecuTec[®] Genu

INNOVATION IN KNEE ORTHOSES

SecuTec[®] Genu Flex SECURELY IN **POSITION WITHOUT** COMPROMISES

69-81 30° 44-52 () 17 cm () 15 cm

BAUERFEIND **MEDICAL AID APP**

Find out more: bauerfeind.de/himi-app





Bauerfeind AG Triebeser Straße 16 07937 Zeulenroda-Triebes Germany

P +49(0)3662866-1000 **F** +49(0)3662866-1999 E info@bauerfeind.com

 \rightarrow BAUERFEIND.COM

SecuTec[®] Genu Flex

SecuTec® Genu Flex

THERAPEUTIC BENEFITS

RELIABLE PROTECTION

... due to the semi-flexible frame that adapts to the leg's shape during movement.

GIVES THE MUSCLES MAXIMUM FREEDOM OF MOVEMENT

> Flexible in the direction of movement and > tension-stable in the direction of therapeutic action.

... with the joint limitable in extension, flexion and immobilization.

 \rightarrow adjustable without tools

THE FIRST FLEXIBLE HARD-FRAME ORTHOSIS

BY BAUERFEIND



"Finally, a 'flexible' orthosis that conforms to the standards of modern medicine as well as the anatomic and individual requirements of the patient."

Christian Teusch Specialist in orthopedics and trauma surgery



The flexible Flex-Links perfectly adapt to movement, remaining securely in position on the leg.

CUSTOMIZED FIT

... through exact adjustment of the circumference by adding or removing Flex-Links.

ightarrow adjustable over the course of treatment





... with extendible hinge legs.

 \rightarrow targeted treatment

HIGH LEVEL OF ACCEPTANCE

... thanks to the flat design and the thin, soft cushioning.

